



PERFIT Tray Adhesive

Convenient and efficient impression accessory

Tray Adhesive for all VPS Materials.


PERFIT Tray Adhesive guarantees excellent adhesion between VPS-based impression materials and metal or plastic trays. The product is remarkably easy to apply and straightforward to use. Impressions can be placed in the tray just 30 minutes after a brief application period of 3 minutes (max) .



VinciSmile PERFIT Tray Adhesive

Indication :

To offer reliable adhesion between impression tray and silicone impression material.



Suitable for A Silicone and C Silicone

Easy to apply and handle




Short setting time after application

Drying Time < 3 minutes

Does not create thickness

Blue colour visual control



Technical Parameters

PERFIT Tray Adhesive

Drying Time	3:00
Bonding Strength (30min)	0.2N/mm ²
Color	Blue
Shelf Life	3 Years

*Data from VinciSmile Lab

Order Information

PERFIT Tray Adhesive

Standard	10ml Bottle
----------	-------------



VINCISMILE GROUP LLC

Add: 20524 Carrey Rd. Walnut, CA 91789, USA

Email: marketing@vincismile.com

Web: www.vincismile.com

Tel: +1 626 283 5808

Shop with us online at store.vincismile.com



Follow @VinciSmile